



# Temperature thermometer

## Introduction

Non-contact Forehead IR Thermometer is specially designed to take the body temperature of a person regardless of room temperature.

Depending on various skin types and thickness, there may be temperature difference.

### 产品参数 / Specification Parameter



Product Name	No contact IR Thermometer
Model No.	LZ600
Display Resolution	0.2°C(0.2°F)
Body Temperature Measurement	34.0-42.9°C / 93-109.4°F
Ambient Temperature Measurement	0-60°C / 32-140°F
Measurement Distance	1 to 15cm(within 5cm is the best)
Auto-shutdown Time	About 7s
Operating Temperature	10-40°C / 50-104°F
Power Supply	2 * 1.5V AAA Batteries(Not Included)
Storage Temperature	0-50°C / 32-122°F
Relative Humidity	≤85%
Item Size	116 * 43 * 152mm (L * W * H)
Item Weight	118g

## Parameters

Product Name: Non Contact IR Thermometer  
Product Model: LZ600  
Product Size: 152 \* 116 \* 43 mm  
Product weight: 118 g  
Body Temperature Range: 34.0 °C-42.9°C, / 93-109°F  
Object Temperature Range: 0 °C -100 °C  
Operating Temperature: 10 °C-40 °C / 50-104°F  
Display Resolution: ±0.2 ° C (0.2°F)  
Measurement Distance: 1-5 cm  
Measurement Time: 0.5 s  
Automatic Power-off function: 15 s  
Memory data: 10 sets Memory  
Relative Humidity: ≤ 85%  
Storage Condition: - 10 °C-40 °C C  
Power Supply: DC 3.0v (2\*AAA batteries)

## Features

Precise non-contact measurements  
User selectable °C or °F  
Selectable Body and surface temp  
Set Alarm Value  
Memorization of the Last 10 measurements  
Automatic Data Hold & Auto power off  
Automatic selection range and  
Display Resolution 0.2°C ( 0.2°F )  
Backlight LCD display.



Product / 产品明细



## SAFETY INFORMATION

This device must only be used for the purposes described in this instruction manual.

This device must only be used in an ambient temperature range between 10 and 40°C.

Do not expose this thermometer to electric shocks.

Do not expose this thermometer to extreme temperature conditions of >50°C or <-20°C.

Do not use the device in relative humidity higher than 85%.

Do not use the device near large electromagnetic fields such as found with cordless or cell phones.

Keep the device away from water and heat, including direct sunlight.

Do not drop or knock the device, and do not use if damaged.

It may affect the accuracy of measurements when the forehead is covered by hair, perspiration, cap or scarf.

Keep the Measuring distance as 1cm-5cm.

Clean the glass with a cotton bud lightly moistened with 75% alcohol.

It may affect the accuracy of measurements when the forehead is covered by perspiration or other factors,

please take the temperature behind the ear lobe.

## Package

Inner packaging: 1pcs / box, box size: 145\*50\*180mm

100pcs/carton, carton size: 750\*500\*370mm weight: 20kg